

# Wound Protector

## Device Description:

- Height-fixed wound protector.
- This device consists of outer ring, inner ring and channel. (See Figure 1)
- Made of silica gel and sterilized.

## Intended Use:

- Expand surgical field in small incision surgery. Avoids incision injury and reduces surgical site infection.

**Table 1: Models and Applicable Incision Size**

Item Code	Size of Incision
HZP1014	03-05 cm
HZP1015	02-04 cm



Figure 1

## Instructions for Use:

- Remove product from sterile package.
- Insert inner ring into the incision. (See Figure 2)
- When removing wound protector, use one hand to grasp the inner ring and carefully remove it from the incision. (See Figure 3)



Figure 2



Figure 3

## Warnings:

- Read all instructions prior to use.
- Do not use if sterile package has been opened or is damaged.
- Re-use or re-sterilization of single-use devices could result in patient morbidity and is improper use of the device.
- Product is for use by a licensed physician only.

## Precautions:

- Product should be smooth and clean.
- Inspect the device for defects prior to use and discard if any defects are found.
- Do not use product if the device expiration date on the label has passed.
- Dispose of the device and packaging in accordance with hospital waste standards and federal regulations.

## Contraindications:

- Not suitable for surgical incision smaller than two centimeters.
- Not suitable for patient allergic to this material.



**BS EN ISO 15223-1**

Symbols used with medical device labels

Symbol	
	Manufacturer
	Authorized representative in the European Community
	Use-by date
	Batch code
	Catalogue number
	Serial number
	Sterile
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not re-sterilize
	Do not use if package is damaged
	Fragile, handle with care
	Keep away from sunlight
	Keep dry
	Temperature limit
	Humidity limitation
	Do not re-use
	Consult instructions for use
	Caution

 For a complete guide of ISO standards visit: [www.iso.org/obp](http://www.iso.org/obp)
**IEC 60601:1**

Symbols for medical electrical equipment

Symbol	Symbol Title
	CLASS II equipment
	TYPE BF APPLIED PART
	General warning sign
	Warning: dangerous voltage
	Refer to instruction manual/booklet <i>*icon may also appear in black, meaning does not change</i>
	WEEE Directive Mark: do not dispose, collect separately <i>IS EN 50419:2006</i>
	IP Classification: Protected against solid foreign objects of $\geq 2,5$ mm. Protection against vertically falling water drops.

Other / Unique Marks to Neurovision Medical Products

Symbol	Symbol Title
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	European Conformity mark
	Not made with natural rubber latex

© 2020 Neurovision Medical Products, Inc. All Rights Reserved.


**Neurovision Medical Products**

 353 Sanjon Road, Ventura, CA, 93001, USA  
[www.neurovisionmedical.com](http://www.neurovisionmedical.com)  
 Tel: 866-815-6999  
 Fax: 877-330-1727

 Jiangsu Haize Medical Scientific  
 Development Co.,Ltd

 Add:No. 99 Furongzhongshan Road  
 Xishan Economic Development Zone  
 Wuzi Jianfu, 214192, China  
 tel: +86-510-82832301 fax: +86-510-82833301

 Shanghai International holding  
 Corp. GmbH (Europe)

 Add: Eiffestrasse 80,20537  
 Hamburg, Germany  
 Tel: +49-40-2513175  
 Fax: +49-40-40-255726  
 Dimid No.: DE/0000040627  
 E-mail: [shholding@hotmail.com](mailto:shholding@hotmail.com)