



RELIABLE UNIVERSAL EMG MONITORING

Unrivaled Precision in Facial Nerve Monitoring

Designed to detect and signal interactions with the facial nerve (CN VII)

- 13 mm needles reach optimal intramuscular location
- 28 G needles designed to reduce skin trauma
- parallel needles are affixed 2.5 mm apart for easy and precise application
- standardized for universal use, inventory management and overall efficiency
- twisted wires and needle insulation protect EMG signals from interference

ORDERING INFORMATION: MASTODON® PAIRED EMG NEEDLE ELECTRODES

Item Code	Description:*
NPDC-5	Mastodon® paired needle electrodes, 2-Channel kit (red & blue), 13 mm, box of 5
NPFC-5	Mastodon® paired needle electrodes, 4-Channel kit (red, orange, blue, purple), 13 mm, box of 5

*each kit includes a STIM return needle and a ground needle