



# SAFE SECURE RETRACTION



Easy Insertion



Dual Ring



Multiple Sizes

## REDUCE SURGICAL SITE INFECTIONS

- Dual-ring, fixed height style protects the tissue & conforms to the incision
- Flexible material inserts easily and maintains its shape for a gentle retraction
- Made with specialized materials, the dual-ring design is atraumatic, durable, and latex-free
- Applicable to various procedures from minimally invasive, small incision (Thyroid, Breast, and Thoracic) to open (GI, C-section, HPB)

**ViClean wound protector** is smooth silica gel that wraps the incision site and protects the surrounding tissue from injury. Dual-ring, fixed height design is atraumatic and expands the surgical site 360° for optimal visualization. **“Bucket” shape self-retracts**, protects the tissue, and keeps the incision clean. Fixed height wound protectors are associated with significantly reducing superficial surgical site infection.<sup>1</sup> ViClean applies to open and minimally invasive surgeries. **More sizes and styles are available**, contact NMP for details.

## ORDERING INFORMATION:

### VICLEAN™ SINGLE-USE WOUND PROTECTOR

Item Code	Description
HZP1014	<b>Wound Protector</b> , Dual-Ring, Fixed Height, Incision Size 2 - 4 cm, Small
HZP1015	<b>Wound Protector</b> , Dual-Ring, Fixed Height, Incision Size 3 - 5 cm, Medium
HZP1016	<b>Wound Protector</b> , Dual-Ring, Fixed Height, Incision Size 4 - 6 cm, Large
HZP1017	<b>Wound Protector</b> , Dual-Ring, Fixed Height, Incision Size 5 - 8 cm, X-Large

1. Zhang, Lisa, Basheer Elsolh, and Sunil V. Patel. "Wound protectors in reducing surgical site infections in lower gastrointestinal surgery: an updated meta-analysis." *Surgical endoscopy* 32.3 (2018): 1111-1122.