

INTRODUCING




PRECISE FLEXIBLE VERSATILE

NERVE IDENTIFICATION & STIMULATION



TAPERED TIP STIMULATION PROBE



- 0.5mm tip diameter tip
- Single-use, monopolar probe
- Ideal for ENT procedures
- Thin, adjustable shaft improves surgical field visibility
- Stimulation for nerve location and tissue mapping

ORDERING INFORMATION: HUMMINGBIRD® SINGLE-USE STIMULATION PROBE

Item Code	Description*:
I-D-TTP90-5	Hummingbird® Tapered Tip Stimulation Probe (Monopolar), 90 mm, Box of 5

**each kit includes subdermal needles, alcohol wipes and NeedleTape®*

SS1-A