



# URETHRAL SPHINCTER ELECTRODE

## Safe. Effective. Needle Free.



- Standard Foley catheter for familiar & reliable function
- Electrodes are flexible & have ZERO profile
- Stimulate & monitor the pudendal/perineal nerve
- Fully integrated EMG capabilities
- No needles or electrode application necessary
- Compatible with standard neurodiagnostic & neuromonitoring systems
- Seamless transition between monitoring & stimulation

### ORDERING INFORMATION:

**USE™ 1-Channel Gold Urethral Catheter Electrode** (each box includes STIM return needle and ground needle)

Item Code	Description
USE8-5	USE™ Gold Urethral Catheter Electrode 8 Fr. <i>(box of 5)</i>
USE10-5	USE™ Gold Urethral Catheter Electrode 10 Fr. <i>(box of 5)</i>
USE12-5	USE™ Gold Urethral Catheter Electrode 12 Fr. <i>(box of 5)</i>
USE14-5	USE™ Gold Urethral Catheter Electrode 14 Fr. <i>(box of 5)</i>
USE16-5	USE™ Gold Urethral Catheter Electrode 16 Fr. <i>(box of 5)</i>
USE18-5	USE™ Gold Urethral Catheter Electrode 18 Fr. <i>(box of 5)</i>

**For more information or to purchase:**




353 Sanjon Rd  
Ventura, CA 93001  
tel: (866) 815-6999  
[www.neurovisionmedical.com](http://www.neurovisionmedical.com)


Product Collaboration by:





# URETHRAL SPHINCTER ELECTRODE



## Needle Free Monitoring of the External Urinary Sphincter



Electrode Detail

-  The USE offers monitoring & stimulation of sacral nerves
-  Wide range of surgical utility:
  - Urology
  - Neuro & spine
  - Rehab & diagnostics applications