



# INTEGRATED SURGICAL PRECISION

new product  
new product  
new product

## Surgical Dissection and Stimulated EMG Combined

- Innovative design ensures a reliable surgical performance and optimizes nerve preservation
- Stimulated EMG enables real-time nerve mapping and continuous monitoring
- Turn-key solution supports informed surgical decisions and efficiency
- Insulated surgical-grade hemostat allows for precise stimulation and tissue dissection
- Universal compatibility with IONM systems for seamless integration
- Value-based pricing ensures affordability without compromising on quality or performance

### ORDERING INFORMATION: SCORPION® SURGICAL INSTRUMENTS

Item Code	Description:*
I-DK-WSF512-5	Scorpion® Monopolar Stimulation Forceps, Fine Curved Tip, Smooth Jaws, 5.5", Single-Use

\*each kit includes a STIM return needle and a ground needle